

www.normanregional.com/lab				
CLIENT LAB SUPPLY FORM				
Order Date:				
Client Name:				
Requested by:				
Address:				
Phone #:				
FAX ALL ORDERS TO 307-1065				
Please order supplies in quantities that will be used before expiration dates. When doing supply inventory, check expiration dates & return supplies that will expire before use.				
Allow 2-4 Business days for delivery				
of in-stock items. (excludes Holidays and inclement weather days)				
To request Supply Forms or check on the status of faxed order, please call 405-307-1150. Supply Forms may be obtained via our website at www.normanregional.com/lab Go to "For Medical Professionals" to print Client Supply Form				
Thank you for your business & the opportunity to assist you in caring for your patients.				
Filled by/Date:				
Delivered by/Date:				
Client/Date:				

Print name:

Backorder	Requested	Order Unit	Supply Name
			Blood Collection
		box of 200	Alcohol Pads
		sets	Blood Culture bottles & Prevantic Swabs
		bag of 250	Cotton balls, medium
		bag of 250	Hubs, Vacutainer, single-use
		box of 100	Needles, syringe, Circle one 21G or 23G
		box of 48	Needles, vacutainer, Circle one 21G or 22G
		each	Syringe - Circle one 3cc, 5cc, or 10cc
		each	Tourniquets
		flat of 100	Tubes, vacutainer, lavender (EDTA, 4.0 ml)
		flat of 100	Tubes, vacutainer, light blue (3.2% Na Citrate, 2.7ml)
		flat of 100	Tubes, vacutainer, gold (gel, 7ml)
		flat of 100	Tubes, vacutainer, mint green (heparin gel, 7ml)
		flat of 100	Tubes, vacutainer, red (plain,7ml)
		each	Tubes, vacutainer, gray (NaF, 7ml)
		each (0)	Tubes, vacutainer, pink (EDTA 6ml)
			Tubes, vacutainer, royal blue (EDTA 6 ml) for lead
		each (2 max)	Tubes, vacutainer, royal blue (Plain)
		o a o b	Non-Blood Specimen Collection
		each	APTIMA® Unisex Swab
		each	APTIMA® Vaginal Swab Specimen Collection Kit (for self-collection)
		each	Eswab Minitip Collection Kit - BD - (green lid)
		each	ESwab <u>Regular</u> Collection Kit - BD - (white lid)
		each	FIT Test Kit (Fecal Immunochemical Kit)
		, ,	Herpes/Viral Culture Media (Universal Transport Media)
		each	MRSA Swab (Copan)
		each	Nasopharyngeal Collection Kit (Flu/RSV/Respiratory Panel)
		, ,	Specimen Containers, Sterile
		, ,	Specipans (Urine/Stool Collection Pans) Stool Container (for Eroch Stool)
		each	Stool Container (for Fresh Stool)
		each each	Stool Kit for C. difficile Toxin - Collection Kit Stool Kit for Cryptosporidium/Giardia (Formalin)
		each	Stool kit for Culture - Collection Kit (Carey-Blair)
		each	Stool Kit for O & P - Collection Kit (LV-PVA/Formalin)
		box	Towelettes, Obstetrical
		each	Throat Culture Swab (BBL)
		each	Urine C&S Transfer Straws
		each	Urine C&S Transfer Tubes (Boric Acid Tubes)
		each	Urine Chemistry Transport Tube - No Additive Vacutainer Tube (white cap)
			Urine Collection & Transport Kit (cup, straw, & tubes)
		each	Vaginosis Panel Kit (Affirm)
	Sr	pecimen Transpo	
		pkg	Specimen Transport Bags
		each	Specimen Transport Bags, STAT (Red)
		each	Specimen Transport Bags, Amber, Light Protected
		Other	
		each	Requisitions - Lab Orders
For Tissue	and Pap	supplies, co	ontact Pathology Consultation Services at 321-0370